

Preparation of 0.1 M Acetate Buffer (pH 5.0)

Aim

To prepare 100 mL of 0.1 M acetate buffer of pH 5.0 using 0.1 M sodium acetate and 0.1 M acetic acid solutions.

Required Solutions

- 0.1 M Acetic acid solution (CH_3COOH)
- 0.1 M Sodium acetate solution (CH_3COONa)
- Distilled water (H_2O)

Apparatus

- 100 mL volumetric flask
- measuring cylinders
- Beaker
- Glass rod
- pH meter

Theory

An acetate buffer consists of a weak acid (acetic acid, CH_3COOH) and its conjugate base (sodium acetate, CH_3COONa). It resists changes in pH when small amounts of acid or base are added. The pH of the buffer is governed by the Henderson–Hasselbalch equation:

$$\text{pH} = \text{p}K_a + \log \left(\frac{[\text{CH}_3\text{COO}^-]}{[\text{CH}_3\text{COOH}]} \right)$$

For acetic acid,

$$\text{p}K_a = 4.76$$

For a buffer of pH 5.0,

$$5.0 = 4.76 + \log \left(\frac{[\text{CH}_3\text{COO}^-]}{[\text{CH}_3\text{COOH}]} \right)$$

Therefore,

$$\frac{[\text{CH}_3\text{COO}^-]}{[\text{CH}_3\text{COOH}]} = 10^{(5.0-4.76)} = 10^{0.24} \approx 1.74$$

Since both stock solutions have the same concentration (0.1 M), the concentration ratio is equal to the volume ratio.

Thus,

$$\frac{V_{\text{CH}_3\text{COONa}}}{V_{\text{CH}_3\text{COOH}}} = 1.74$$

For a total volume of 100 mL,

$$V_{\text{CH}_3\text{COONa}} = \frac{1.74}{1 + 1.74} \times 100 = 63.5 \text{ mL}$$

$$V_{\text{CH}_3\text{COOH}} = \frac{1}{1 + 1.74} \times 100 = 36.5 \text{ mL}$$

Hence, mixing 63.5 mL of 0.1 M sodium acetate with 36.5 mL of 0.1 M acetic acid produces a 0.1 M acetate buffer of pH 5.0.

Procedure

1. Clean a 100 mL volumetric flask and rinse it with distilled water.
2. Measure 63.5 mL of 0.1 M sodium acetate (CH_3COONa) using measuring cylinder and add into the volumetric flask.
3. Add 36.5 mL of 0.1 M acetic acid (CH_3COOH) into above solution.
4. Mix the solution thoroughly with constant stirring.
5. Measure the pH using a calibrated pH meter.
6. If necessary, adjust the pH:
 - Add a few drops of 0.1 M acetic acid to decrease the pH.
 - Add a few drops of 0.1 M sodium acetate to increase the pH.
7. If total volume is slightly less than 100 mL after pH adjustment, make up the volume to mark with distilled water.

Observation

The prepared solution was clear and colourless. The measured pH was approximately 5.0.

Result

A 100 mL sample of 0.1 M acetate buffer of pH 5.0 was successfully prepared by mixing 63.5 mL of 0.1 M sodium acetate solution with 36.5 mL of 0.1 M acetic acid solution.

Precautions

1. Wash the Apparatus carefully.
2. pH meter must be calibrated before use.
3. Handle Acid solution with care.